

# Video Nystagmus - Brille

## Video Nystagmus Spectacles




Innovative Video – Nystagmusbrille mit zwei hochwertigen Videochips für eine optimale Dunkeluntersuchung des Nystagmus.

The Video Nystagmography Spectacles allows an optimum examination in the dark with high resolution IR-cameras.

Mit dieser Video – Nystagmusbrille wird die Beobachtung des Nystagmus ideal ermöglicht. Sie zeichnet sich durch eine hohe Qualität, **hohe Lebensdauer** und ein optimales Preis-/Leistungsverhältnis aus.

Durch die **anpassungsfähige, sehr bequeme und lichtundurchlässige**

**Brille** wird der Patient minimal belastet und seine Augen in **vollkommene Dunkelheit** gehüllt. In der Brille sorgen **zwei hochwertige Infrarot-Kameras** für eine perfekte Wiedergabe auf dem Monitor. Der **leistungsstarke und wartungsfreie Akku** liefert die Spannung für das gesamte System und kann durch das **mitgelieferte Ladegerät** einfach und schnell wieder aufgeladen werden.



Video images of the eyes are obtained without direct contact by this system. This system offers high quality, the best possible price/performance relation. Because the goggles are **flexible, convenient and impervious to light** the patient are burdened minimally and the eyes are covered in **darkness**

**completely.**

**Two high quality infrared cameras** ensure a good image on the screen.

The accumulator is **powerful** and **maintenance free** and provide all parts of the system by voltage. The accumulator can be charged quick and easily by supplied charger.

### Technische Daten:

Artikelnummer: 08-415 (Set, komplett)  
 Stromversorgung: Akku-Pack (08-415-30)  
 Batterielaufzeit: ca. 2 h

#### Kameramodule

Auflösung: 352 x 288  
 Öffnungswinkel: 56° x 42°  
 Brennweite: 4,9  
 Blende: 2,8

#### Akku-Pack

Blei-Gel

Spannung: 12 V  
 Kapazität: 1300 mAh

#### Ladegerät

Ladespannung: 13,8 V

#### Monitor

Display: TFT LCD, 7"  
 Betriebsspannung: 12 V  
 Leistungsaufnahme: 8 W

### Technical data:

Item code: 08-415 (Set, complete)  
 Power source: Accumulator (08-415-30)  
 Battery run time: Approximately 2 h

#### Cameras

Resolution: 352 x 288  
 Field of View: 56° x 42°  
 Focal length: 4.9  
 Aperture: 2.8

#### Accumulator

Lead-gel  
 Voltage: 12 V  
 Capacity: 1300 mAh

#### Charger

Voltage: 13.8 V

#### Monitor

Display method: TFT LCD, 7"  
 Operating voltage: 12 V  
 Power consumption: 8 W