

## Platin-Iridium-Brenner

Die wechselbaren Platin-Iridium-Brenner sind alle autoklavierbar und für die Anwendung mit unserem leistungsstarken Kaustiktransformator optimal geeignet.

Durch die Verwendung von Platin ist gewährleistet, daß die Brenner über die gesamte Lebensdauer nicht oxidieren (keine Schwärzung!). Der Zusatz von Iridium bewirkt Härte und lange Lebensdauer.

Das breite Spektrum an Brennerformen und -längen eröffnet dem Anwender vielfältige Nutzungsmöglichkeiten.

## Tips of Platinum-Iridium

The interchangeable tips of Platinum-Iridium are autoclavable and fit optimal to our high performance Cautery Transformer.

Using platinum is the reason for non oxidation of the tips over lifetime (no blackening!). Owing to the addition of iridium, these cauteries are distinguished by their extreme hardness and durability.

The wide spectrum of different shapes and lengths of tips enables to use them for many different purposes.

Im vielfältigen praktischen Einsatz haben sich die folgenden Brennerlängen als sinnvoll erwiesen:

Auge: 3 cm

Zahn/Haut: 5 cm

Mund/Haut: 10 cm

Ohr/Nase: 15 cm

Kehlkopf/Nasenrachen: 20 cm

Over many years of practical use it was found that following lengths of Tips are most helpful:

Auge: 3 cm

Dental / Skin: 5 cm

Mouth / Skin: 10 cm

Ear / Nose: 15 cm

Larynx / Nasopharynx: 20 cm

Form/Shape	Länge/Length				
	3 cm	5 cm	10 cm	15 cm	20cm
<b>A</b>	09-161	09-201	09-221	09-241	09-261
<b>B</b>	09-162	09-202	09-222	09-242	09-262
<b>C</b>	09-163	09-203	09-223	09-243	09-263
<b>D</b>	09-164	09-204	09-224	09-244	09-264
<b>E</b>	09-165	09-205	09-225	09-245	09-265
<b>F</b>	09-166	09-206	09-226	09-246	09-266
<b>G</b>	09-167	09-207	09-227	09-247	09-267
<b>H</b>	09-168	09-208	09-228	09-248	09-268
<b>I</b>	09-169	09-209	09-229	09-249	09-269
<b>K</b>	09-170	09-210	09-230	09-250	09-270
<b>L</b>	09-171	09-211	09-231	09-251	09-271
<b>M</b>	09-172	09-212	09-232	09-252	09-272
<b>N</b>	09-173	09-213	09-233	09-253	09-273
<b>O</b>	09-174	09-214	09-234	09-254	09-274
<b>P</b>	09-175	09-215	09-235	09-255	09-275

