

## Ersatzlampen passend für Martin-Leuchten

### Replacement lamps fitting Martin OT-lights

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting modell	Spannung Voltage	Watt Wattage	Sockel Base	Bild Pic.
99-61230	86-903-00		24 V	55 W	G 6.35	2
99-61250	86-905-01		24 V	120 W	G 6.35	2
99-61270	86-907-01		24 V	150 W	G 6.35	2
99-61210	86-910-00		24 V	250 W	GY 6.35	2
99-61252	88-905-24		24 V	100 W	G 6.35	2
99-61262	88-906-24		24 V	150 W	G 6.35	2
99-61207	88-907-01		24 V	150 W	G 6.35	2
99-61272	88-907-24		24 V	250 W	G 6.35	2
99-61200	88-910-02		24 V	100/100 W	E27	-
99-61220	88-912-02		24 V	75/75 W	E27	-
99-61240	88-914-02		24 V	50/50 W	BA22D	-
99-61260	88-916-02		24 V	50 W	BA20D	3
99-61214	88-916-04		24 V	25 W	BA20D	3
99-61298	98-004-7007		15 V	150 W	Litzen	-

## Ersatzlampen passend für Berchtold-Leuchten und Projektoren


### Replacement lamps fitting Berchtold OT-lights and projectors

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base	Bild Pic.
99-612K1	100052579	Kaltlichtprojektor / Cold light source	15 V	150 W	Litzen	
99-612C2	CZ 904-24	C-200 / C-300 / D-300	24 V	55 W	G 6.35	2
99-612C4	CZ 905-22	C-450 / C-452	22.8 V	110 W	G 6.35	2
99-61252	CZ 905-24	C-450alt / B-480 / B-680 / B-950	24 V	100 W	G 6.35	2
99-61262	CZ 906-24	C-570	22.8 V	150 W	G 6.35	2
99-612C9	CZ 907-22	C-950G / C-952, D-650	22.8 V	250 W	G 6.35	2
99-61210	CZ 907-24	C 950	24 V	250 W	GY 6.35	2
99-612C6	CZ 908-22		22.8 V	150 W	G 6.35	2
99-612C5	CZ 908-24	C-571 / C-572	24 V	150 W	G 6.35	2
99-61260	CZ 916-02		24 V	50 W	BA20D	3
99-61214	CZ 916-04		24 V	25 W	BA20D	3
99-612D2	CZ 940-12	D-200	12 V	35 W	G 6.35	2
99-612D3	CZ 955-12	D-300	12 V	50 W	G 6.35	2
99-612K2	CZ 960-12		12 V	20 W	GU 5.3	

## Ersatzlampen passend für Amsco-Leuchten

### Replacement lamp fitting Amsco OT-lights


Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting modell	Spannung Voltage	Watt Wattage	Sockel Base	Bild Pic.
99-63508	P-093902-608		22 V	50 W		-
99-63501	P-129249-001		33 V	235 W	GY9.5	-
99-63528	P-129362-228		22 V	220 W	GY9.5	-
99-63534	P-129359-934		20 V	200 W	G17	-



1



2



3