

Ersatzlampen passend für Leitz-Mikroskope
Replacement lamps fitting Leitz microscopes

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base
99-74717	37-717		24 V	150 W	G 6.35
99-74723	37-723		24 V	250 W	G 6.35
99-74520	050-520			200 W	SFc10-4
99-74535	050-535		6 V	10 W	BA15S
99-74550	050-550		6 V	10 W	G 4
99-74610	050-610			200 W	SFc10-4
99-74615	050-615			200 W	SFc10-4
99-74008	500-008		20 V	150 W	SFcX12-4
99-74010	500-010	Panphot, Metallux ND, Orthoplan	12 V	60 W	BA20D
99-74011	500-011	Monla-Leuchte, Laborlux 2	6 V	5 A	E 14
99-74012	500-012	Mikrospektrograph, Dilatometer	6 V	15 W	BA15D
99-74016	500-013			200 W	SFc10-4
99-74031	500-031		17 V	450 W	SK19/36
99-74073	500-073	HM-Lux3, HM-Lux, SM-Lux	6 V	5 W	BA9S
99-74096	500-096	SM-Lux, TL-Leuchte von HM-Lux	6 V	10 W	BA15s
99-74137	500-137	LH 100Z/102Z/103Z/105Z	39-45 V	50 W	Sfa6-2
99-74138	500-138	LH 100Z/102Z/103Z/105Z	20 V	100 W	Sfa7,5-2
99-74139	500-139		14 V	75 W	Sfa7,5-2
99-74153	500-153	LH 250/252	61 V	200 W	SFc10-4
99-74177	500-177	Dialux	6 V	15 W	P25d
99-74182	500-182	LH 50, Ortholux 2, Epivert	12 V	50 W	G 6.35
99-74245	500-245	SM-Lux, Laborlux, Ergolux, u.a.	6 V	20 W	G4
99-74910	500-910	Promar, Prado-Univers., Pradovit	24 V	250 W	G 6.35
99-74974	500-974	LH 100/102/103/104/105	12 V	100 W	G 6.35
99-74620	505-620		61 V	200 W	SFc10-4
99-74025	610-025		6 V	4,5 A	E14/16
99-74071	620-071		6 V	5 A	Ba20d
99-74334		OP-Mikroskop Leica M680	12 V	50 W	Spezial
99-74606			21 V	150 W	G5.3
99-S022	11063	Leica Bilirubinmeter, model 10311C			

Ersatzlampen passend für Wild-Mikroskope
Replacement lamps fitting Wild microscopes

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base
99-74324	166324		6 V	15 W	BA 15d
99-74354	166354/000113			200 W	SFc10-4
99-74355	166355		20 V	150 W	SFcX12-4
99-76760	177160		6 V	20 W	P21D
99-74049	193049		12 V	100 W	GY 6.35
99-74643	384643		12 V	50 W	Spezial
99-74053	997053		15 V	150 W	GZ 6.35
99-74084	997084		24 V	250 W	G 6.35
99-74334	10445334		12 V	50 W	Spezial

Ersatzlampen passend für Zeiss-Jena Mikroskope
Replacement lamps fitting Zeiss-Jena microscopes

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base
99-71147			6 V	25 W	PY16-1.25
99-71133	214797	Laboval 4	6 V	20 W	PY16-1.25
99-71612	LWT-F4		12 V	50 W	Ba20d
99-71251	LWT-P2		6 V	15 W	M16
99-71272	LWT-P5		6 V	15 W	M16
99-71273	LWT-P5T		6 V	15 W	M16
99-Y1058	17161		6 V	5 W	Ba15s