

Ersatzlampen passend für Olympus-Geräte
Replacement lamps fitting Olympus devices

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base
99-75342	035-342		120 V	30 W	Ba15d
99-75120	1120 / 1121		6 V	20 W	G 4
99-75815	300815		15 V	150 W	GZ6.35
99-75624	624		6 V	10 W	G 4
99-S007	8-B-152	CH2	120 V	30 W	Ba15d
99-75B90	8-B-190			200 W	Sfc10-4
99-75B92	8-B-192		20 V	100 W	Sfa9-2
99-75B93	8-B-193		23 V	50 W	Sfa8-2
99-75C11	8-C-101	LSO, LSE, MTX	6 V	5 A	Ba15d
99-75S30	8-C-102	BHA, BHB, IMT	6-8 V	5 A	special
99-75S15	8-C-103	BMMJ, STM, BHC, BHM, VM-LSG	6 V	2,5 A	special
99-75C01	8-C-401	A-LSH-3, BH-RFL, B-LSH-3, PME-LSM	12 V	100 W	G 6.35
99-75C02	8-C-402	BHM-SLH-2, BHM-MLSH	12 V	50 W	G 6.35
99-75C03	8-C-403		6 V	10 W	G4
99-75C04	8-C-404		21 V	150 W	GX 5.3
99-75C05	8-C-405		6 V	20 W	G4
99-75C06	8-C-406		6 V	10 W	G4
99-75G01	8-G-101		6 V	5 W	E10
99-75102	8-G-102		6 V	10 W	E10
99-75M02	8-M-102		6 V	15 W	BA15D
99-75S01	8-S-401		15 V	150 W	GZ 6.35
99-75S06	8-S-006		20 V	150 W	SfcX12-4
99-75134	PG 55134		12 V	20 W	Litzen
99-S080	SB 0493		220 V	20 W	Ba15d
99-75500	ZM 036500		11 V	50 W	15d/Flansch

Ersatzlampen passend für Nikon-Mikroskope
Replacement lamps fitting Nikon microscopes

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting Model	Spannung Voltage	Watt Wattage	Sockel Base
99-72249	70249		6 V	20 W	P15D
99-72410	70410				
99-72818	71818		10 V	2 A	P28D
99-72826	71826		24 V	150 W	G 6.35
99-72168	72168		15 V	150 W	GZ 6.35
99-72223	76223		21 V	150 W	GX 5.3
99-72451	76451		20 V	150 W	SfCX12-4
99-72513	76513		12 V	100 W	G 6.35
99-72516	76516		12 V	50 W	G 6.35
99-72906	76906		6 V	15 W	Ba20d
99-72907	76907				
99-72098	77098		6 V	0.04 A	Midget
99-72903	77903		6 V	15 W	P25D
99-72904	77904		6 V	30 W	BA15D
99-72910	77910		6 V	15 W	P25D
99-72153	78153		17 V	450 W	SK19/36
99-72564	78564		39...45 V	50 W	Sfa6-2
99-72580	78580			200 W	SFf10-4
99-72590	78590			200 W	SFf10-4
99-72099	79099		6 V	20 W	G 4
99-72727	79727		12 V	5 W	Bi-Pin
99-72542	MXA-20-542		240 V	30 W	
99-72316	MXA-20-316		6 V	15 W	Ba15d
99-720SE		SE	220 V	20 W	Ba15d
99-Y5761		Eclipse E400	6 V	30 W	G 4
99-S126	500WS	Xenon Blitzlampe NF-505			